



# IVM-C3 (Confocal v. 3)

## The New All-in-One Intravital Imaging Platform



### Tractable, Fast and Gentle

IVM-C3 is a highly integrated intravital microscopy for *in vivo* imaging with an enormously increased detection efficiency, optical resolution, and contrast of the image compared to conventional fluorescence microscopy. Equipped with a 4-wavelength laser and four high-sensitivity confocal detectors, IVM-C3 allows multi-dimensional views of living cells and tissues in 3D or 4D up to four different colors. It is the optimal system for simultaneously observing various dynamic multi-cellular behaviors in small live animal models.



### Key Features

- Multi-color Simultaneous Imaging (4 channels, 4 different colors)
- Fully Integrated *in vivo* Maintenance Unit / Animal Stage (e.g., Monitoring & Homeostatic Regulation of Animal Vitality)
- Ultra High-Speed Imaging (max. 100 fps - 512x512 pixels)
- 4D Animal Motion Compensation (XYZ & Time)


# IVM-C3 (Confocal v. 3)

## The New All-in-One Intravital Imaging Platform

### SPECIFICATIONS

<b>Laser</b>	<b>Confocal Laser Unit</b>	<ul style="list-style-type: none"> <li>• 405 nm (20mW), 488 nm (20mW), 561 nm (20mW), 640 nm (20mW)</li> </ul>
<b>Fluorescence Detector</b>	<b>Confocal Detector</b>	<ul style="list-style-type: none"> <li>• Wavelength : 185-900 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>• 4 Ultra-broadband high SNR PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> <li>• 25-2,000 <math>\mu\text{m}</math> variable pinhole</li> </ul>
	<b>Variable Emission Filter (Optional)</b>	<ul style="list-style-type: none"> <li>• 2 or 6 emission filters can be mounted on each of four detectors</li> </ul>
<b>Scan Head</b>	<b>Scanner</b>	<ul style="list-style-type: none"> <li>• Polygonal mirror (Fast axis scanning, Max. 66 kHz)</li> <li>• Galvano scanner (Slow axis scanning, Max. 200 <math>\mu\text{s}/\text{step}</math>)</li> </ul>
<b>Imaging Head</b>	<b>Objectives</b>	<ul style="list-style-type: none"> <li>• Max. 5 objectives are mountable on S/W controlled motorized turret (1X – 100X)</li> <li>• Compatible for commercial objectives</li> </ul>
<b>Image</b>	<b>FOV</b>	<ul style="list-style-type: none"> <li>• 100 x 100 <math>\mu\text{m}^2</math> - 10 x 10 <math>\text{mm}^2</math></li> </ul>
	<b>Pixel Resolution</b>	<ul style="list-style-type: none"> <li>• Max. 2,048 x 2,048 pixels</li> </ul>
	<b>Imaging Speed</b>	<ul style="list-style-type: none"> <li>• Standard : 30 fps @ 512 x 512 pixels</li> <li>• (Optional) High Speed : 60 fps @ 512 x 512 pixels</li> <li>• (Optional) Ultra HighSpeed : 100 fps @ 512 x 512 pixels</li> </ul>
<b>Animal / Sample Stage</b>	<b>3D Stage</b>	<ul style="list-style-type: none"> <li>• Travel Range : 50,000 x 50,000 x 75,000 <math>\mu\text{m}</math> (XYZ)</li> <li>• Micromanipulation (Max. 0.2 <math>\mu\text{m}</math> resolution)</li> <li>• 3-axis independent control with Jog Dial &amp; IVM Engine software</li> </ul>
	<b>Specimen Holder</b>	<ul style="list-style-type: none"> <li>• Flexible-design universal <i>in vivo</i> / <i>ex vivo</i> / <i>in vitro</i> specimen holder can be mounted</li> <li>• (Optional) Homeothermic warming system, Holders for window chamber</li> </ul>
<b>Animal Motion Compensation</b>	<b>4D <i>In vivo</i> Imaging Motion Compensation</b>	<ul style="list-style-type: none"> <li>• XY motion compensation : Averaged image acquisition with motion artifact compensation</li> <li>• Z motion compensation : Image-based sample Z position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>• T motion compensation : Image-based image XY position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>• Combination of above three compensation for 4D <i>in vivo</i> motion compensation</li> <li>• Controllable by IVM Engine software</li> </ul>
<b>Additional <i>In vivo</i> Modules</b>	<b>Live Animal Maintenance Unit</b>	<ul style="list-style-type: none"> <li>• Body Temp. Monitoring &amp; Feedback Heater Control, including tablet PC</li> <li>• 4CH Rectal Probe, Body Plate Heater, Thermometer Sensor &amp; Cover Glass Heater</li> </ul>
	<b><i>In vivo</i> Imaging Chamber SET</b>	<ul style="list-style-type: none"> <li>• Standard Dorsal Skinfold Chamber SET</li> <li>• Lung Imaging Chamber SET</li> <li>• Cranial Window SET</li> <li>• Abdominal Imaging Window SET</li> <li>• Pancreas Imaging Window SET</li> <li>• Mammary Imaging Window SET</li> <li>• Heart Imaging Window SET</li> </ul>
	<b>Inhalation Anesthesia System</b>	<ul style="list-style-type: none"> <li>• Whole Rodent Animal Inhalation Anesthesia System</li> <li>• Anesthesia Mask and Connections</li> </ul>
<b>Engine &amp; Studio Software</b>	<b>Image Display</b>	<ul style="list-style-type: none"> <li>• Independent 4 single channel display (RGBA channel)</li> <li>• Overlay channel display (Selection among RGBA channel)</li> </ul>
	<b><i>In vivo</i> Imaging Mode</b>	<ul style="list-style-type: none"> <li>• Mosaic imaging (XY), Z-stackimaging (Z), Time-lapse imaging (T)</li> <li>• Time-lapse imaging at Multi-position (T - M),</li> <li>• Time-lapse &amp; Z-stack imaging (TZ),</li> <li>• Time-lapse &amp; Z-stack imaging at Multi-position (TZ - M)</li> </ul>


### New All-in-One IVM Series Size Information



**IVM Technology, Inc.**  
**Webpage** [www.ivimtech.com](http://www.ivimtech.com)  
**Contact** [information@ivimtech.com](mailto:information@ivimtech.com)  
**TEL** +82-2-431-7450  
**FAX** +82-2-3400-0450



# IVM-M3 (Two-Photon v. 3)

## The New All-in-One Intravital Imaging Platform



### Deep Tissue Imaging, High-resolution, Tunable Laser

IVM-M3 is one of the renowned members of The New All-in-One IVM Series. It has the flexibility feature of the traditional converted microscope and the high-resolution imaging ability of second-harmonic generation microscopy. It is equipped with a fully-automated tunable fs-pulse NIR laser system. IVM-M3 is the optimal system for users who need to conduct deeper-tissue imaging using less-scattering NIR wavelength. Full control functionality of the fs-laser system is integrated with the two-photon imaging software for user convenience with various automation algorithms.



### Key Features

- Deep tissue imaging with a tunable long-wavelength NIR fs-laser system
- Fully Integrated *in vivo* Maintenance Unit / Animal Stage (e.g., Monitoring & Homeostatic Regulation of Animal Vitality)
- Ultra High-Speed Imaging (max. 100 fps - 512x512 pixels)
- 4D Animal Motion Compensation (XYZ & Time)
- Label-free, non-linear second and third harmonic generation ability


# IVM-M3 (Two-Photon v. 3)

## The New All-in-One Intravital Imaging Platform

### SPECIFICATIONS

<b>Laser</b>	<b>Tunable Two-Photon Laser Unit</b>	<ul style="list-style-type: none"> <li>• Ti : Sapphire laser</li> <li>• Wavelength : 690-1050 nm, Pulse width &lt; 75 fs, Rep. rate : 80 MHz</li> <li>• Avg. power &gt; 2.5 W, Dispersion compensation : 0 to -49,000 fs<sup>2</sup></li> </ul>
<b>Fluorescence Detector</b>	<b>Two-Photon Detector</b>	<ul style="list-style-type: none"> <li>• Wavelength : 185 - 760 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>• 4 High quantum efficiency PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> </ul>
	<b>Variable Emission Filter (Optional)</b>	<ul style="list-style-type: none"> <li>• 2 or 6 emission filters can be mounted on each of four detectors</li> </ul>
<b>Scan Head</b>	<b>Scanner</b>	<ul style="list-style-type: none"> <li>• Polygonal mirror (Fast axis scanning, Max. 66 kHz)</li> <li>• Galvano scanner (Slow axis scanning, Max. 200 <math>\mu</math>s / step)</li> </ul>
<b>Imaging Head</b>	<b>Objectives</b>	<ul style="list-style-type: none"> <li>• Max. 5 objectives are mountable on S/W controlled motorized turret (1X - 100X)</li> <li>• Compatible for commercial objectives</li> </ul>
<b>Image</b>	<b>FOV</b>	• 100 x 100 $\mu$ m <sup>2</sup> - 10 x 10 mm <sup>2</sup>
	<b>Pixel Resolution</b>	• Max. 2,048 x 2,048 pixels
	<b>Imaging Speed</b>	<ul style="list-style-type: none"> <li>• Standard : 30 fps @ 512 x 512 pixels</li> <li>• (Optional) High Speed : 60 fps @ 512 x 512 pixels</li> <li>• (Optional) Ultra High Speed : 100 fps @ 512 x 512 pixels</li> </ul>
<b>Animal / Sample Stage</b>	<b>3D Stage</b>	<ul style="list-style-type: none"> <li>• Travel Range : 50,000 x 50,000 x 75,000 <math>\mu</math>m (XYZ)</li> <li>• Micromanipulation (Max. 0.2 <math>\mu</math>m resolution)</li> <li>• 3-axis independent control with Jog Dial &amp; IVM Engine software</li> </ul>
	<b>Specimen Holder</b>	<ul style="list-style-type: none"> <li>• Flexible-design universal <i>in vivo</i> / <i>ex vivo</i> / <i>in vitro</i> specimen holder can be mounted</li> <li>• (Optional) Homeothermic warming system, Holders for window chamber</li> </ul>
<b>Animal Motion Compensation</b>	<b>4D <i>In vivo</i> Imaging Motion Compensation</b>	<ul style="list-style-type: none"> <li>• XY motion compensation : Averaged image acquisition with motion artifact compensation</li> <li>• Z motion compensation : Image-based sample Z position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>• T motion compensation : Image-based image XY position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>• Combination of above three compensation for 4D <i>in vivo</i> motion compensation</li> <li>• Controllable by IVM Engine software</li> </ul>
<b>Additional <i>In vivo</i> Modules</b>	<b>Live Animal Maintenance Unit</b>	<ul style="list-style-type: none"> <li>• Body Temp. Monitoring &amp; Feedback Heater Control, including tablet PC</li> <li>• 4CH Rectal Probe, Body Plate Heater, Thermometer Sensor &amp; Cover Glass Heater</li> </ul>
	<b><i>In vivo</i> Imaging Chamber SET</b>	<ul style="list-style-type: none"> <li>• Standard Dorsal Skinfold Chamber SET</li> <li>• Lung Imaging Chamber SET</li> <li>• Cranial Window SET</li> <li>• Abdominal Imaging Window SET</li> <li>• Pancreas Imaging Window SET</li> <li>• Mammary Imaging Window SET</li> <li>• Heart Imaging Window SET</li> </ul>
	<b>Inhalation Anesthesia System</b>	<ul style="list-style-type: none"> <li>• Whole Rodent Animal Inhalation Anesthesia System</li> <li>• Anesthesia Mask and Connections</li> </ul>
<b>Engine &amp; Studio Software</b>	<b>Image Display</b>	<ul style="list-style-type: none"> <li>• Independent 4 single channel display (RGBA channel)</li> <li>• Overlay channel display (Selection among RGBA channel)</li> </ul>
	<b><i>In vivo</i> Imaging Mode</b>	<ul style="list-style-type: none"> <li>• Mosaic imaging (XY), Z-stack imaging (Z), Time-lapse imaging (T)</li> <li>• Time-lapse imaging at Multi-position (T- M),</li> <li>• Time-lapse &amp; Z-stack imaging (TZ),</li> <li>• Time-lapse &amp; Z-stack imaging at Multi-position (TZ- M)</li> </ul>

### New All-in-One IVM Series Size Information




IVM Technology, Inc. Webpage [www.ivimtech.com](http://www.ivimtech.com) | Contact [information@ivimtech.com](mailto:information@ivimtech.com)  
TEL +82-2-431-7450 | FAX +82-2-3400-0450




# IVM-CM3 (Confocal and Two-Photon v. 3)

## The New All-in-One Intravital Imaging Platform



### High Contrast and Resolution, Dual-mode and Tunable Laser

IVM-CM3 is a highly integrated All-in-One intravital microscopy. It can focus on the desired wavelength with its tunable Two-Photon laser unit for wavelengths as low as 690 nm, higher up to 1050 nm, or in between. IVM-CM3 combines the advantages of both Confocal and Two-Photon microscopy providing endless possibilities for three-dimensional imaging of living cells near the skin or deep into the tumor in small animals.



### Key Features

- Deep tissue imaging with a tunable long-wavelength NIR fs-laser system
- One-click automated transition between Confocal and Two-Photon imaging modes
- Fully Integrated *in vivo* Maintenance Unit / Animal Stage (e.g., Monitoring & Homeostatic Regulation of Animal Vitality)
- Ultra High-Speed Imaging (max. 100 fps - 512x512 pixels)
- 4D Animal Motion Compensation (XYZ & Time)


# IVM-CM3 (Confocal and Two-Photon v. 3)

## The New All-in-One Intravital Imaging Platform

### SPECIFICATIONS

Laser	Confocal Laser Unit	<ul style="list-style-type: none"> <li>405 nm (20mW), 488 nm (20mW), 561 nm (20mW), 640 nm (20mW)</li> </ul>
	Tunable Two-Photon Laser Unit	<ul style="list-style-type: none"> <li>Ti : Sapphire laser</li> <li>Wavelength : 690-1050 nm, Pulse width &lt; 75 fs, Rep. rate : 80 MHz</li> <li>Avg. power &gt; 2.5 W, Dispersion compensation : 0 to -49,000 fs<sup>2</sup></li> </ul>
Fluorescence Detector	Confocal Detector	<ul style="list-style-type: none"> <li>Wavelength : 185 - 900 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>4 Ultra-broadband high SNR PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> <li>25-2,000 <math>\mu</math>m variable pinhole</li> </ul>
	Two-Photon Detector	<ul style="list-style-type: none"> <li>Wavelength : 185 - 760 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>4 High quantum efficiency PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> </ul>
	Variable Emission Filter (Optional)	<ul style="list-style-type: none"> <li>2 or 6 emission filters can be mounted on each of four detectors</li> </ul>
Scan Head	Scanner	<ul style="list-style-type: none"> <li>Polygonal mirror (Fast axis scanning, Max. 66 kHz)</li> <li>Galvano scanner (Slow axis scanning, Max. 200 <math>\mu</math>s /step)</li> </ul>
Imaging Head	Objectives	<ul style="list-style-type: none"> <li>Max. 5 objectives are mountable on S/W controlled motorized turret (1X – 100X)</li> <li>Compatible for commercial objectives</li> </ul>
Image	FOV	100 x 100 $\mu$ m <sup>2</sup> - 10 x 10 mm <sup>2</sup>
	Pixel Resolution	Max. 2,048 x 2,048 pixels
	Imaging Speed	<ul style="list-style-type: none"> <li>Standard : 30 fps @ 512 x 512 pixels</li> <li>(Optional) High Speed : 60 fps @ 512 x 512 pixels</li> <li>(Optional) Ultra High Speed : 100 fps @ 512 x 512 pixels</li> </ul>
Animal / Sample Stage	3D Stage	<ul style="list-style-type: none"> <li>Travel Range : 50,000 x 50,000 x 75,000 <math>\mu</math>m (XYZ)</li> <li>Micromanipulation (Max. 0.2 <math>\mu</math>m resolution)</li> <li>3-axis independent control with Jog Dial &amp; IVM Engine software</li> </ul>
	Specimen Holder	<ul style="list-style-type: none"> <li>Flexible-design universal <i>in vivo</i> / <i>ex vivo</i> / <i>in vitro</i> specimen holder can be mounted</li> <li>(Optional) Homeothermic warming system, Holders for window chamber</li> </ul>
Animal Motion Compensation	4D <i>In vivo</i> Imaging Motion Compensation	<ul style="list-style-type: none"> <li>XY motion compensation : Averaged image acquisition with motion artifact compensation</li> <li>Z motion compensation : Image-based sample Z position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>T motion compensation : Image-based image XY position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>Combination of above three compensation for 4D <i>in vivo</i> motion compensation</li> <li>Controllable by IVM Engine software</li> </ul>
Additional <i>In vivo</i> Modules	Live Animal Maintenance Unit	<ul style="list-style-type: none"> <li>Body Temp. Monitoring &amp; Feedback Heater Control, including tablet PC</li> <li>4CH Rectal Probe, Body Plate Heater, Thermometer Sensor &amp; Cover Glass Heater</li> </ul>
	<i>In vivo</i> Imaging Chamber SET	<ul style="list-style-type: none"> <li>Standard Dorsal Skinfold Chamber SET</li> <li>Lung Imaging Chamber SET</li> <li>Cranial Window SET</li> <li>Abdominal Imaging Window SET</li> <li>Pancreas Imaging Window SET</li> <li>Mammary Imaging Window SET</li> <li>Heart Imaging Window SET</li> </ul>
	Inhalation Anesthesia System	<ul style="list-style-type: none"> <li>Whole Rodent Animal Inhalation Anesthesia System</li> <li>Anesthesia Mask and Connections</li> </ul>
Engine & Studio Software	Image Display	<ul style="list-style-type: none"> <li>Independent 4 single channel display (RGBA channel)</li> <li>Overlay channel display (Selection among RGBA channel)</li> </ul>
	<i>In vivo</i> Imaging Mode	<ul style="list-style-type: none"> <li>Mosaic imaging (XY), Z-stack imaging (Z), Time-lapse imaging (T)</li> <li>Time-lapse imaging at Multi-position (T- M),</li> <li>Time-lapse &amp; Z-stack imaging (TZ),</li> <li>Time-lapse &amp; Z-stack imaging at Multi-position (TZ- M)</li> </ul>

### New All-in-One IVM Series Size Information



IVM Technology, Inc. Webpage [www.ivimtech.com](http://www.ivimtech.com) | Contact [information@ivimtech.com](mailto:information@ivimtech.com)  
TEL +82-2-431-7450 | FAX +82-2-3400-0450



# IVM-MS3 (Two-Photon Smart v. 3)

## The New All-in-One Intravital Imaging Platform



### Compact, Cost-saving, Hands-free

IVM-MS3 is the smart version of IVM-M3, an All-in-One Two-Photon Intravital Microscopy optimized for *in vivo* imaging. It integrates a compact, high-stability and maintenance-free fs-pulse laser unit into a single box. IVM-MS3 is perfectly capable of imaging deep tissues within a wavelength fixed at 920nm, which makes it an excellent resource for researchers with a specific target but limited resources and budget.



### Key Features

- Simple and hands-free turn-key operation of 920 nm NIR fs-laser for deep tissue imaging
- Label-free, non-linear second and third harmonic generation ability
- Fully Integrated *in vivo* Maintenance Unit / Animal Stage (e.g., Monitoring & Homeostatic Regulation of Animal Vitality)
- Ultra High-Speed Imaging (max. 100 fps – 512x512 pixels)
- 4D Animal Motion Compensation (XYZ & Time)




# IVM-MS3 (Two-Photon Smart v. 3)

## The New All-in-One Intravital Imaging Platform

### SPECIFICATIONS

<b>Laser</b>	<b>Compact Two-Photon Laser Unit</b>	<ul style="list-style-type: none"> <li>Air cooled fs-fiber laser system with built-in power control</li> <li>Wavelength : 920 nm, Pulse width &lt;150 fs, Rep. rate : 80 MHz</li> <li>Avg. power &gt;0.8 W, Dispersion compensation : 0 to -22,000 fs<sup>2</sup> ○</li> </ul>
<b>Fluorescence Detector</b>	<b>Two-Photon Detector</b>	<ul style="list-style-type: none"> <li>Wavelength : 185 - 760 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>4 High quantum efficiency PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> </ul>
	<b>Variable Emission Filter (Optional)</b>	<ul style="list-style-type: none"> <li>2 or 6 emission filters can be mounted on each of four detectors</li> </ul>
<b>Scan Head</b>	<b>Scanner</b>	<ul style="list-style-type: none"> <li>Polygonal mirror (Fast axis scanning, Max. 66 kHz)</li> <li>Galvano scanner (Slow axis scanning, Max. 200 μs /step)</li> </ul>
<b>Imaging Head</b>	<b>Objectives</b>	<ul style="list-style-type: none"> <li>Max. 5 objectives are mountable on S/W controlled motorized turret (1X – 100X)</li> <li>Compatible for commercial objectives</li> </ul>
<b>Image</b>	<b>FOV</b>	<ul style="list-style-type: none"> <li>100 x 100 μm<sup>2</sup> - 10 x 10 mm<sup>2</sup></li> </ul>
	<b>Pixel Resolution</b>	<ul style="list-style-type: none"> <li>Max. 2,048 x 2,048 pixels</li> </ul>
	<b>Imaging Speed</b>	<ul style="list-style-type: none"> <li>Standard : 30 fps @ 512 x 512 pixels</li> <li>(Optional) High Speed : 60 fps @ 512 x 512 pixels</li> <li>(Optional) Ultra High Speed : 100 fps @ 512 x 512 pixels</li> </ul>
<b>Animal / Sample Stage</b>	<b>3D Stage</b>	<ul style="list-style-type: none"> <li>Travel Range : 50,000 x 50,000 x 75,000 μm (XYZ)</li> <li>Micromanipulation (Max. 0.2 μm resolution)</li> <li>3-axis independent control with Jog Dial &amp; IVM Engine software</li> </ul>
	<b>Specimen Holder</b>	<ul style="list-style-type: none"> <li>Flexible-design universal <i>in vivo</i> / <i>ex vivo</i> / <i>in vitro</i> specimen holder can be mounted</li> <li>(Optional) Homeothermic warming system, Holders for window chamber</li> </ul>
<b>Animal Motion Compensation</b>	<b>4D <i>In vivo</i> Imaging Motion Compensation</b>	<ul style="list-style-type: none"> <li>XY motion compensation : Averaged image acquisition with motion artifact compensation</li> <li>Z motion compensation : Image-based sample Z position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>T motion compensation : Image-based image XY position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>Combination of above three compensation for 4D <i>in vivo</i> motion compensation</li> <li>Controllable by IVM Engine software</li> </ul>
<b>Additional <i>In vivo</i> Modules</b>	<b>Live Animal Maintenance Unit</b>	<ul style="list-style-type: none"> <li>Body Temp. Monitoring &amp; Feedback Heater Control, including tablet PC</li> <li>4CH Rectal Probe, Body Plate Heater, Thermometer Sensor &amp; Cover Glass Heater</li> </ul>
	<b><i>In vivo</i> Imaging Chamber SET</b>	<ul style="list-style-type: none"> <li>Standard Dorsal Skinfold Chamber SET</li> <li>Lung Imaging Chamber SET</li> <li>Cranial Window SET</li> <li>Abdominal Imaging Window SET</li> <li>Pancreas Imaging Window SET</li> <li>Mammary Imaging Window SET</li> <li>Heart Imaging Window SET</li> </ul>
	<b>Inhalation Anesthesia System</b>	<ul style="list-style-type: none"> <li>Whole Rodent Animal Inhalation Anesthesia System</li> <li>Anesthesia Mask and Connections</li> </ul>
<b>Engine &amp; Studio Software</b>	<b>Image Display</b>	<ul style="list-style-type: none"> <li>Independent 4 single channel display (RGBA channel)</li> <li>Overlay channel display (Selection among RGBA channel)</li> </ul>
	<b><i>In vivo</i> Imaging Mode</b>	<ul style="list-style-type: none"> <li>Mosaic imaging (XY), Z-stack imaging (Z), Time-lapse imaging (T)</li> <li>Time-lapse imaging at Multi-position (T- M),</li> <li>Time-lapse &amp; Z-stack imaging (TZ),</li> <li>Time-lapse &amp; Z-stack imaging at Multi-position (TZ- M)</li> </ul>

### New All-in-One IVM Series Size Information




**IVM Technology, Inc.**  
**Webpage** [www.ivimtech.com](http://www.ivimtech.com)  
**Contact** [information@ivimtech.com](mailto:information@ivimtech.com)  
**TEL** +82-2-431-7450  
**FAX** +82-2-3400-0450




# IVM-CMS3 (Confocal and Two-Photon Smart v. 3)

## The New All-in-One Intravital Imaging Platform



### Cost-Effective, Straightforward, Dual-mode

IVM-CMS3 is the world's most compact and affordable dual-mode intravital confocal and two-photon microscope, providing versatile functionality in a single box. Having the Confocal laser units of IVM-C3 and the compact Two-Photon laser unit of IVM-MS3 with a one-switch mode changing feature, IVM-CMS3 provides comfortable multi-purpose use for intravital functional imaging and saves from unnecessary space and high costs.



### Key Features

- Simple and hands-free turn-key operation of 920 nm NIR fs-laser for deep tissue imaging
- One-click automated transition between Confocal and Two-Photon imaging modes
- Fully Integrated *in vivo* Maintenance Unit / Animal Stage (e.g., Monitoring & Homeostatic Regulation of Animal Vitality)
- Ultra High-Speed Imaging (max. 100 fps - 512x512 pixels)
- 4D Animal Motion Compensation (XYZ & Time)
- Label-free, non-linear second and third harmonic generation ability


# IVM-CMS3 (Confocal and Two-Photon Smart v. 3)

## The New All-in-One Intravital Imaging Platform

### SPECIFICATIONS

Laser	Confocal Laser Unit	<ul style="list-style-type: none"> <li>405 nm (20mW), 488 nm (20mW), 561 nm (20mW), 640 nm (20mW)</li> </ul>
	Compact Two-Photon Laser Unit	<ul style="list-style-type: none"> <li>Air cooled fs-fiber laser system with built-in power control</li> <li>Wavelength : 920 nm, Pulse width &lt;150 fs, Rep. rate : 80 MHz</li> <li>Avg. power &gt;0.8 W, Dispersion compensation : 0 to -22,000 fs<sup>2</sup></li> </ul>
Fluorescence Detector	Confocal Detector	<ul style="list-style-type: none"> <li>Wavelength : 185 - 900 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>4 Ultra-broadband high SNR PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> <li>25-2,000 <math>\mu</math>m variable pinhole</li> </ul>
	Two-Photon Detector	<ul style="list-style-type: none"> <li>Wavelength : 185 - 760 nm (DAPI, CFP, GFP, YFP, RFP, Cy5, Cy5.5, etc.)</li> <li>4 High quantum efficiency PMTs (UV to Near IR, Ultra High Sensitivity, Low Dark Current)</li> </ul>
	Variable Emission Filter (Optional)	<ul style="list-style-type: none"> <li>2 or 6 emission filters can be mounted on each of four detectors</li> </ul>
Scan Head	Scanner	<ul style="list-style-type: none"> <li>Polygonal mirror (Fast axis scanning, Max. 66 kHz)</li> <li>Galvano scanner (Slow axis scanning, Max. 200 <math>\mu</math>s /step)</li> </ul>
Imaging Head	Objectives	<ul style="list-style-type: none"> <li>Max. 5 objectives are mountable on S/W controlled motorized turret (1X - 100X)</li> <li>Compatible for commercial objectives</li> </ul>
Image	FOV	<ul style="list-style-type: none"> <li>100 x 100 <math>\mu</math>m<sup>2</sup> - 10 x 10 mm<sup>2</sup></li> </ul>
	Pixel Resolution	<ul style="list-style-type: none"> <li>Max. 2,048 x 2,048 pixels</li> </ul>
	Imaging Speed	<ul style="list-style-type: none"> <li>Standard : 30 fps @ 512 x 512 pixels</li> <li>(Optional) High Speed : 60 fps @ 512 x 512 pixels</li> <li>(Optional) Ultra High Speed : 100 fps @ 512 x 512 pixels</li> </ul>
Animal / Sample Stage	3D Stage	<ul style="list-style-type: none"> <li>Travel Range : 50,000 x 50,000 x 75,000 <math>\mu</math>m (XYZ)</li> <li>Micromanipulation (Max. 0.2 <math>\mu</math>m resolution)</li> <li>3-axis independent control with Jog Dial &amp; IVM Engine software</li> </ul>
	Specimen Holder	<ul style="list-style-type: none"> <li>Flexible-design universal <i>in vivo</i> / <i>ex vivo</i> / <i>in vitro</i> specimen holder can be mounted</li> <li>(Optional) Homeothermic warming system, Holders for window chamber</li> </ul>
Animal Motion Compensation	4D <i>In vivo</i> Imaging Motion Compensation	<ul style="list-style-type: none"> <li>XY motion compensation : Averaged image acquisition with motion artifact compensation</li> <li>Z motion compensation : Image-based sample Z position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>T motion compensation : Image-based image XY position adjustment for long-term intravital microscopic imaging &amp; sample tracking (Feedback-loop automatic stage control)</li> <li>Combination of above three compensation for 4D <i>in vivo</i> motion compensation</li> <li>Controllable by IVM Engine software</li> </ul>
Additional <i>In vivo</i> Modules	Live Animal Maintenance Unit	<ul style="list-style-type: none"> <li>Body Temp. Monitoring &amp; Feedback Heater Control, including tablet PC</li> <li>4CH Rectal Probe, Body Plate Heater, Thermometer Sensor &amp; Cover Glass Heater</li> </ul>
	<i>In vivo</i> Imaging Chamber SET	<ul style="list-style-type: none"> <li>Standard Dorsal Skinfold Chamber SET</li> <li>Lung Imaging Chamber SET</li> <li>Cranial Window SET</li> <li>Abdominal Imaging Window SET</li> <li>Pancreas Imaging Window SET</li> <li>Mammary Imaging Window SET</li> <li>Heart Imaging Window SET</li> </ul>
	Inhalation Anesthesia System	<ul style="list-style-type: none"> <li>Whole Rodent Animal Inhalation Anesthesia System</li> <li>Anesthesia Mask and Connections</li> </ul>
Engine & Studio Software	Image Display	<ul style="list-style-type: none"> <li>Independent 4 single channel display (RGBA channel)</li> <li>Overlay channel display (Selection among RGBA channel)</li> </ul>
	<i>In vivo</i> Imaging Mode	<ul style="list-style-type: none"> <li>Mosaic imaging (XY), Z-stack imaging (Z), Time-lapse imaging (T)</li> <li>Time-lapse imaging at Multi-position (T- M),</li> <li>Time-lapse &amp; Z-stack imaging (TZ),</li> <li>Time-lapse &amp; Z-stack imaging at Multi-position (TZ- M)</li> </ul>

### New All-in-One IVM Series Size Information



**IVM Technology, Inc.**  
**Webpage** [www.ivimtech.com](http://www.ivimtech.com)  
**Contact** [information@ivimtech.com](mailto:information@ivimtech.com)  
**TEL** +82-2-431-7450  
**FAX** +82-2-3400-0450